## PRODUCT DESCRIPTION

**HARD BITE REGISTRATION MATERIAL** is a standard bite registration material based on vinylpolysiloxanes suitable for all types of bite registrations.

**BITE REGISTRATION MATERIAL** is a bite registration material with a mousse-like consistency based on vinylpolysiloxanes suitable for all types of bite registrations.

**CLEAR BITE REGISTRATION MATERIAL** is a transparent material based on vinylpolysiloxanes suitable for all types of bite registrations that require transparent matrices. These transparent matrices can be used for the veneering technique.

TECHNICAL DATA	HARD BITE REGISTRATION MATERIAL	BITE REGISTRATION MATERIAL			CLEAR BITE REGISTRATION MATERIAL
ISO 4823		Type B : bite registration materials			
Indications	<ul style="list-style-type: none"> <li>Accurate occlusal registrations</li> <li>Key material for needle point registration</li> <li>Production of small model segments</li> <li>For applications that require hard bite registrations</li> </ul>	<ul style="list-style-type: none"> <li>Accurate occlusal registrations</li> <li>Key material for needle point registration</li> <li>Production of small model segments</li> </ul>			<ul style="list-style-type: none"> <li>Accurate occlusal registrations</li> <li>Key material for needle point registration</li> <li>Production of small model segments</li> <li>Silicone key for composite restorations</li> </ul>
		Super Fast Set	Fast Set	Regular Set	
REF	FAST SET BR-9010	SUPER FAST SET BR-9003 BR-9005	FAST SET BR-9004 BR-9002	REGULAR SET BR-9001	TRANSPARENT BR-9200
*Total working time (minutes)	0:40	0:30	0:40	1:30	1:00
Time in mouth (minutes)	1:00	0:45	1:00	3:00	2:00
*Total setting time (minutes)	1:40	1:15	1:40	4:30	3:00
Linear dimensional change (%)	<0.1	<0.1			<0.3
Shore A hardness	93±2	85±2			70±2

## ATTENTION

Application is at 23°C ± 2°C / 73°F ± 4°F, 50 ± 5% relative humidity. Higher temperatures reduce working and setting times and lower temperatures prolong them.

Intended Use: Dental Impression material

Use environment: Dental practice

Patient target group: People treated in the course of a dental procedure with no limitations in age or gender.

## MOUNTING THE CARTRIDGE

Insert cartridge onto dispensing gun following the manufacturer's instructions and remove cap. Extrude impression material until it exits both chambers at the same rate and then mount the appropriate mixing tip. An intra-oral tip may be placed on the mixing tip.

**APPLICATION - BITE REGISTRATIONS:** Apply bite registration material according to the desired technique (occlusal or buccal). Extrude material first onto areas where teeth have been prepared keeping the end of the mixing tip immersed in the material to avoid the formation of air bubbles and ask the patient to bite. Carefully remove the set bite registration after the recommended time. Remove all residues of the material from the mouth. Trim bite registration to ensure correct seating on model. After use keep mixing tip on cartridge as it serves as a cap. Intraoral tips and mixing tips are for single use only.

**APPLICATION - TRANSPARENT MATRIX:** Apply CLEAR BITE REGISTRATION MATERIAL onto appropriate impression tray keeping the end of the mixing tip immersed in the material to avoid the formation of air bubbles and place in mouth. Carefully remove the set matrix after the recommended time. Rinse with water and blow dry with air. Prepare teeth according to the desired technique and following the manufacturer's instructions of use of the composite material. Apply composite onto the transparent matrix and position in mouth. Light cure through the matrix following the corresponding instructions of use. After use keep mixing tip on cartridge as it serves as a cap. Intraoral tips and mixing tips are for single use only.

**DISINFECTING:** Once impression is removed from mouth, rinse under lukewarm water. The impression can then be disinfected with appropriate disinfecting solutions or sprays, such as 2% glutaraldehyde. Take care to follow the manufacturer's instructions for use. For further information regarding disinfecting impressions: JADA 1991; vol. 122; IS(8); p110. If disinfecting with hydrogen peroxide, thoroughly rinse with lukewarm water to avoid the formation of bubbles.

**PRECAUTIONS:** • Contact with latex gloves may impair the setting reaction. Vinyl gloves are recommended • Eugenol or certain haemostatic materials may impair the setting reaction • Do not combine with condensation curing silicones, polyether or polysulfide • Avoid spillage of material on clothing as it cannot be removed chemically • Avoid contact of material with eyes. In case of contact rinse eyes with plenty of water and consult a physician immediately • Seek immediate medical attention in case of ingestion • Do not use after the use-by date

**CONTRAINDICATIONS:** Application of the above material is contraindicated:

- If a patient is known to be hypersensitive or allergic to any of the ingredients
- If a patient has loose teeth, since they can be further loosened or extracted by taking an impression

**STORAGE:** Store product between 10-25°C (50-77°F). Store impression at room temperature (max. 25°C/77°F).

**LOT NUMBER AND EXPIRY DATE:** The above are printed on each package and should be quoted in all correspondence which requires identification of the product.

**DISPOSAL:** Dispose of the product according to local regulations as applicable

**USERS / FOR PROFESSIONAL USE ONLY:** These materials have been developed for use exclusively by dental professionals and must be handled according to the instructions of use. The recommendations included in the instructions correspond to the actual state of knowledge regarding dental techniques. This does not relieve the user from determining the suitability of the product for its intended purpose. The manufacturer cannot be held responsible for damages caused by other uses or incorrect handling.

Caution: Federal USA law restricts this device to sale by or on the order of a dental professional.

**REPORTING OBLIGATIONS:** Should a serious incident in relation to this product occur, it should be reported to the manufacturer and the competent authority.

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