6575 Butler Creek Road, Missoula, MT 59808
(800) 551-5172 • Fax (406) 549-7452 • www.am-eagle.com

UK Responsible Person
Astek Innovations Ltd
95 Atlantic Street, Altrincham
WA14 5DH United Kingdom
Tel: +44 (0)161 942 3900

EC REP

Young Microbrush Ireland Ltd.
Clogherane, Dunganarvan
Co. Waterford,
X35 VE02, Ireland

UK CA CH REP

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

USA

MADE IN USA

NON-STERILE

MD

MEDICAL DEVICE

DOUBLE GRACEY KITS	
AEDGKXPX	AEGEK2XPX
AEDGMKXPX	AEGEK6XPX
AEGEKXPX	

EN - Dental Instrument - Use for intended purpose only. Inspect, clean and sterilize prior to each use. Do not heat above 300° F/148° C.
ES - Instrumento dental - Usar únicamente para la finalidad indicada. Examine, limpie y esterilize antes de cada uso. No calentar sobre 300° F/148° C.
DE - Zahnmedizinisches Instrument - Nur für den vorgesehenen Zweck verwenden. Vor jedem Gebrauch inspizieren, reinigen und sterilisieren. Nicht über 300° F/148° C erhitzen.
FR - Instrument dentaire - Destiné à un usage spécifique. Contrôler, nettoyer et stériliser avant chaque utilisation. Ne pas chauffer à plus de 300° F/148° C.
IT - Strumento dentale - Solo per uso specificato, Ispezionare, pulire e sterilizzare prima di ciascun uso, Non riscaldare a piu di 300° F/148° C.

EN INSTRUMENT CARE

Each AEI instrument has been designed and manufactured with the greatest of care to fulfill specific functional criteria. Incorrect handling or misuse reduces the service life of these precision instruments. Designed use is thus an integral part of proper instrument care.

Sterilization cannot be a substitute for cleaning! Stubborn impurities and debris should be removed with a soft brush (never with steel wool, drill brushes or abrasive items). The instrument should be rinsed in distilled water (not tap water). Some dental surgeries incorporate ultrasonic units or thermo-disinfectors within their cleaning systems, where extreme care is also required to ensure that the cleaning agents and chemicals being used do not damage the instrument subsurface. These methods clean but do not sterilize!

AEI recommends that instruments be disinfected using solutions which are phenol-free.

Instruments should be cleaned in a non-corrosive, neutral cleaning agent with minimal foaming.

AEI recommends using cassettes for sterilizing any of your instruments. Cassettes offer the protection your instruments need to avoid damage by contact with other instruments, especially when using ultrasonic cleaning systems. They also offer proper flow of air or water for more efficient cleaning and sterilization.

Do not batch stainless steel, aluminum, brass or copper instruments together during the cleaning or sterilization processes. If batched together, a potential for electrolysis reactions between dissimilar metals will exist, which can produce etching and corrosion on the instrument surfaces.

Careful drying of instruments during cleaning and sterilizing processes is extremely important. Any remaining water or condensation can cause potential rust or corrosion on the instrument surface and into the substrate.

With proper care, the longevity and performance of all your professional dental instruments can be extended. We hope that the above has been helpful to you in achieving this goal.

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