



# Biozyme-LT™

Low Temperature Enzymatic Instrument Presoak & Cleaner

## BIOTROL™ BIOZYME-LT™ ENZYMATIC INSTRUMENT PRE-SOAK AND CLEANER

Biotrol™ Biozyme-LT™ is a concentrated dual-enzyme detergent designed to dissolve blood and tissue on contaminated instruments within 2 to 5 minutes. It works effectively in both room temperature and warm water, reducing excessive instrument handling and manual scrubbing. Ideal for use as a pre-soak or in ultrasonic cleaners.

### Key Features and Benefits:

**Fast-Acting Formula:** Dissolves organic debris in 2 to 5 minutes.

**Dual-Temperature Efficiency:** Effective in both warm and room temperature water.

**Low Foaming:** Suitable for ultrasonic use without excessive suds.

**Built-In Corrosion Inhibitor:** Protects instrument surfaces during the cleaning process.

**Minty Fresh Scent:** Enhances user experience with a pleasant fragrance.

**Usage Instructions:** Dilute by mixing 1 oz. of concentrate per 1 gallon of water (1:128 ratio).

Product Specifications	
Item Number	BLT900CS (4ct)
Size	30 oz. (32 oz.) bottle
Dilution Ratio	1:128 (1 oz. per gallon of water)
Fragrance	Mint

### Contact:

Visit [younginnovations.com/biotrol](http://younginnovations.com/biotrol) or contact your local Dealer representative for more details and ordering assistance.

## ORDER TODAY!

Available through authorized Biotrol™ dealers.

