



Ultrasonic Enzymatic Tablets

A Powerful and Simple Cleaning Solution

DEFEND® Ultrasonic Enzymatic Tablets are a convenient and effective solution for cleaning dental instruments. These easy-to-use tablets are formulated with fast-acting protease enzymes, which assist in the removal of blood, bone, tissue, food particles and other soils that may be resistant to chemical detergents. The biodegradable, non-chlorine, non-corrosive formulation ensures safety for both instruments and users.



Key Features and Benefits:

Rapid Dissolution: Tablets dissolve in water in under a minute, streamlining the cleaning process.

Effective Cleaning: Enzymatic action targets and breaks down organic debris, enhancing cleaning efficacy.

Versatile Use: Suitable for use as an ultrasonic cleaning solution, instrument pre-soak, or evacuation system cleaner.

User-Friendly: Eliminates the need to handle and store heavy liquid concentrates.

Pleasant Fragrance: Features a neutral pH and a pleasant scent for improved user experience.



ORDER TODAY!

Available through your authorized DEFEND dealers

Product Specifications	
Item Number	UT-1000
Quantity	64 tablets/box
Dilution Ratio	Two tablets per gallon of water
Formulation	Biodegradable, non-chlorine, non-corrosive with fast-acting protease enzymes
pH Level	Neutral